

Nasal Shell

The Nasal Shell is an excellent reconstruction option for correcting nasal deformities. The implant mimics the shape of the nasal bones and upper lateral cartilage. Nasal Shell inserts can be used in dorsal areas where additional augmentation is required. An identical sterile implant is included for use as a template. Sutures may be placed in any desired location.

Cat#	Description	A	B	C							
4368	Regular	38mm	21mm	17mm							
	Insert - Small (included)	30mm	9mm	4mm							
	Insert - Large (included)	38mm	9mm	2.5mm							
4369	Large	40mm	18mm	Insert - Small (included)	32mm	9mm	4mm	Insert - Large (included)	41mm	9mm	3.1mm
	Insert - Small (included)	32mm	9mm	4mm							
	Insert - Large (included)	41mm	9mm	3.1mm							

